

Harloff MR Conditional 7 Drawer Anesthesia Cart, Specialty Pkg, Key Lock | MR7K-MAN



Stock Photo - Options May Vary

FACTORY NEW

MR-Conditional Seven Drawer Anesthesia Cart Key Lock Specialty Package MR7K-MAN. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware.

The cart is Latex-Free and has a Hard-baked powder-coated finish. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites.

Tested by an independent testing authority to 3 Tesla levels. Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional."

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- 13 drawer colors available (see below)

SPECIFICATIONS

- Assembled Cart Dimensions: 68.5" H x 32" W x 25" D
- Cart net weight: 137 pounds
- Approx. Shipping weight: 157 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)