

Harloff MR-Cond. Narrow 7 Dr. Anesthesia Cart w/Standard Pkg & Key Lock | MRN7K



Stock Photo - Options May Vary

FACTORY NEW

The Harloff MRN7K MR-Conditional Narrow Seven Drawer Anesthesia Cart Key Lock Standard Package. Construction: Non-magnetic or “weakly magnetic” materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware.

The cart is Latex Free and has a Hard-baked powder-coated finish. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites.

Tested by an independent testing authority to 3 Tesla levels. Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional.”

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI

SPECIFICATIONS

- Assembled Cart Dimensions: 446” H x 20” W x 23.5” D
- Cart net weight: 95 pounds
- Approx. Shipping weight: 110 pounds

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com