

# Harloff MR-Conditional 6 Drawer Anesthesia Cart w/Standard Pkg. & Key Lock | MR6K



Stock Photo - Options May Vary

## FACTORY NEW

Harloff MR6K MRI-Conditional Six Drawer Anesthesia Cart Key Lock Standard Package. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Additionally, the Cart is Latex-Free and has a Hard-baked powder-coated finish.

Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites. Tested by an independent testing authority to 3 Tesla levels. Furthermore.

meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional."

## FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full-extension drawer slides

## SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 30.5" W x 23.75" D
- Cart net weight: 126 pounds
- Approx. Shipping weight: 146 pounds

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)