

Harloff MR-Conditional 7 Drawer Anesthesia Cart w/Standard Package & Key Lock | MR7K



Stock Photo - Options May Vary

FACTORY NEW

Harloff MR7K MR-Conditional Seven Drawer Anesthesia Cart Key Lock Standard Package. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. The cart is Latex Free and has a Hard-baked powder-coated finish.

Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites. Tested by an independent testing authority to 3 Tesla levels. Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional."

Standard Features Include Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full-extension drawer slides

SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 32" W x 25" D
- Drawer Configuration: (1) - 4", (4) - 3.25", (1) - 6.5", and (1) - 9.75"
- Cart net weight: 137 pounds
- Approx. Shipping weight: 157 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*