

Harloff MR-Conditional Narrow 6 Drawer Anesthesia Cart w/Std Pkg & Key Lock | MRN6K



Stock Photo - Options May Vary

FACTORY NEW

Harloff 6 Drawer Anesthesia Cart w/Std Pkg & Key Lock, model number MRN6K. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. The cart is Latex-Free and has a Hard-baked powder-coated finish.

Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites and tested by an independent testing authority to 3 Tesla levels. Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional."

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI identification sticker applied

SPECIFICATIONS

- Assembled Cart Dimensions: 46" H x 20" W x 23.5" D
- Cart net weight: 95 pounds
- Approx. shipping weight: 110 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)