

Harloff MR-Conditional Narrow 7 Drawer Anesthesia Cart w/Specialty Pkg | MRN7K-MAN



Stock Photo - Options May Vary

FACTORY NEW

Harloff Narrow 7 Drawer Anesthesia Cart w/Specialty Pkg: MRN7K-MAN. This 7 Drawer Harloff MR-Conditional Cart with Specialty Package has been designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites.

An independent testing authority has tested it to 3 Tesla levels. Additional features include compliance with ASTM International, Designation: F2503-05, which defines the cart as 'MR-Conditional.' The cart's narrow design makes it ideal for MRI suites with limited floor space.

The cabinet is painted gloss white and features an MRI identification sticker for easy recognition. Users can choose from 13 available drawer colors, offering customization options to suit various needs and preferences.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Narrow design for use in MRI suite with limited floor space
- Cabinet painted gloss white with MRI identification sticker applied

SPECIFICATIONS

- MR Conditional
- Cabinet: painted white w/MRI ID sticker applied
- Finish: Hard-Baked Powder Coat
- Drawer Slides: Low-ferrous, ball-bearing, full-extension drawer slides
- Casters: Low-ferrous 5" sealed ball-bearing casters; full swivel - one with brake, one-directional

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)