

Harloff MR-Conditional 6 Drawer ER Cart, Standard Pkg. w/Breakaway Lock | MR6B



Stock Photo - Options May Vary

FACTORY NEW

Harloff 6 Drawer ER Cart, Standard Pkg w/Breakaway Lock, model MR6B. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites and tested by an independent testing authority to 3 Tesla levels.

Furthermore, this unit Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional".

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- 13 drawer colors available (see below)

SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 30.5" W x 23.75" D
- Internal Drawer Dimensions: 23.9" W x 17" D
- Drawer Configuration: (1) - 4", (2) - 3.25", (2) - 6.5", and (1) - 9.75"
- Cart Weight: Net Weight: 126 pounds
- Approx. Shipping Weight: 146 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)