

Harloff MR-Conditional 6 Drawer ER Cart w/Breakaway Lock Specialty Pkg | MR6B-EMG



Stock Photo - Options May Vary

FACTORY NEW

The Harloff 6 Drawer ER Cart with Breakaway Lock Specialty Package (model MR6B-EMG) is constructed from non-magnetic and weakly magnetic materials including aluminum, stainless steel, and plastic. It meets ASTM International Designation F2503-05 for MR-Conditional use and has been independently tested to 3 Tesla MRI field strength, making it safe for use in MRI suites.

The cart features low-ferrous ball-bearing full-extension drawer slides and 5" sealed ball-bearing casters with full swivel capability. A stainless steel breakaway locking system provides security, while a full wrap-around vinyl bumper and protective corner bumpers protect equipment and facility surfaces. The cabinet is painted gloss white with an MRI identification sticker and includes a replaceable plastic top with integrated push handles.

The Emergency Specialty Package includes stainless steel IV poles, a cardiac board with mounting brackets, an oxygen tank holder for D/E cylinders, four aluminum utility hooks, and two full drawer trays with adjustable dividers. The cart is latex-free with a hard-baked powder coat finish and comes fully assembled with pre-threaded accessory mounting holes. It includes a twelve-year warranty.

FEATURES

- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Low-ferrous, ball-bearing, full-extension drawer slides
- Low-ferrous 5" sealed ball-bearing casters; full swivel - 1 with brake, and one directional
- Stainless steel IV poles (MR-IV2)
- Twelve-year warranty

SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 32" W x 25" D
- Internal Drawer Dimensions: 23.9" W x 17" D
- Drawer Configuration: (1) - 4", (2) - 3.25", (2) - 6.5", and (1) - 9.75"
- IV Pole: Extends to 68.75" H
- Cart Weight: Net Weight: 138 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com