

Harloff MR-Conditional 7 Drawer Emergency Cart Standard Pkg | MR7B



Stock Photo - Options May Vary

FACTORY NEW

Harloff 7 Drawer Emergency Cart Standard Package, model MR7B. MR-Conditional Six Drawer Emergency Cart Breakaway Lock Specialty Package MR6B-EMG.

Construction: Non-magnetic or “weakly magnetic” materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites, and tested by an independent testing authority to 3 Tesla levels.

Additionally, it Meets the American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional. Standard Features Include: Cabinet painted gloss white with MRI identification sticker applied, Low-ferrous, ball-bearing, and full extension drawer slides.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full-extension drawer slides

SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 32" W x 25" D
- Internal Drawer Dimensions: 23.9" W x 17" D
- Drawer Configuration: (1) - 4", (4) - 3.25", (1) - 6.5", and (1) - 9.75"
- Cart net weight: 138 lbs
- Approx. Shipping weight: 158 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)