

# Harloff MR-Conditional 7 Drawer Emergency Cart Specialty Pkg | MR7B- EMG



Stock Photo - Options May Vary

## FACTORY NEW

Harloff 7 Drawer Emergency Cart Specialty Package, model MR7B-EMG. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Furthermore, it is designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites and tested by the independent testing authority to 3 Tesla level.

## FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Cabinet painted gloss white with MRI identification sticker applied
- 13 drawer colors available

## SPECIFICATIONS

- Assembled Cart Dimensions: 45.75" H x 32" W x 25" D
- Drawer Configuration: (1) - 4", (4) - 3.25", (1) - 6.5", and (1) - 9.75"
- Internal Drawer Dimensions: 23.9" W x 17" D
- IV pole: extends to 68.75" H
- Cart net weight: 138 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)