

Harloff MR-Conditional Narrow 6 Drawer ER Cart w/Breakaway Lock Standard Pkg | MRN6B



Stock Photo - Options May Vary

FACTORY NEW

Harloff Narrow 6 Drawer ER Cart w/Breakaway Lock Standard Pkg, model MRN6B. Construction: Non-magnetic or “weakly magnetic” materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites. Tested by an independent testing authority to 3 Tesla levels. Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional. Narrow design for use in MRI suites where floor space is an issue.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI identification sticker applied

SPECIFICATIONS

- Assembled Cart Dimensions: 46” H x 20” W x 23.5” D
- Internal Drawer Dimensions: 13.75” W x 17” D
- Drawer Configuration: (1) - 4”, (2) - 3.25”, (2) - 6.5”, and (1) - 9.75”
- Cart Weight: Net Weight: 95 pounds
- Approx. Shipping Weight: 110 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*