

Harloff MR-Conditional Narrow 6 Dr ER Cart w/Breakaway Lock Specialty Pkg | MRN6B-EMG



Stock Photo - Options May Vary

FACTORY NEW

Harloff MRN6B-EMG MR-Conditional Narrow 6 Drawer ER Cart w/Breakaway Lock Specialty Pkg. 20 Years of Innovation in Crash Cart Design and Construction has resulted in the Universal Line.

Its modern appearance incorporates industry-leading standard features with flexible configuration options. Completely equipped for demanding Emergency situations with a complete accessory package.

Furthermore, this cart has a durable powder-coated paint finish over aluminum construction.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI identification sticker applied

SPECIFICATIONS

- Assembled Cart Dimensions: 46" H x 20" W x 23.5" D
- Drawer Dimensions: 13.75 W x 17" D
- Drawer Configuration: (1) - 4", (2) - 3.25", (2) - 6.5", and (1) - 9.75"
- IV Pole: Extends to 68.75" H
- Cart Weight: Net Weight: 113 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)