

Harloff MR-Conditional Narrow 7 Drawer ER Cart w/Breakaway Lock Standard Pkg | MRN7B



Stock Photo - Options May Vary

FACTORY NEW

Harloff 7 Drawer ER Cart w/Breakaway Lock Standard Pkg, model MRN7B. Construction: Non-magnetic or "weakly magnetic" materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware. Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites and tested by an independent testing authority to 3 Tesla levels.

Furthermore, the unit meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional. Standard features include a narrow design for use in MRI suites with limited floor space. The cabinet is painted gloss white with an MRI identification sticker.

Choose from 13 available drawer colors.

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI identification sticker applied

SPECIFICATIONS

- Assembled Cart Dimensions: 46" H x 20" W x 23.5" D
- Internal Drawer Dimensions: 13.75" W x 17" D
- Drawer Configuration: (1) - 4", (4) - 3.25", (1) - 6.5", (1) - 9.75"
- Cart net weight: 95 pounds
- Approx. Shipping weight: 110 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)