

Harloff MR-Conditional Narrow 7 Drawer ER Cart w/Breakaway Lock Specialty Pkg | MRN7B-EMG



Stock Photo - Options May Vary

FACTORY NEW

Harloff Narrow 7 Drawer ER Cart w/Breakaway Lock Specialty Pkg, model MRN7B-EMG. Construction: Non-magnetic or “weakly magnetic” materials, including aluminum, stainless steel, plastic, and aluminum and brass attachment hardware.

Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites, and tested by an independent testing authority to 3 Tesla levels. Furthermore, this unit meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of “MR-Conditional.”

FEATURES

- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Tested by an independent testing authority to 3 Tesla level
- Narrow design for use in MRI suites where floor space is an issue
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full-extension drawer slides

SPECIFICATIONS

- Assembled Cart Dimensions: 46” H x 20” W x 23.5” D
- IV Pole: Extends to 68.75” H
- Internal Drawer Dimensions: 13.75” W x 17” D
- Drawer Configuration: (1) - 4”, (4) - 3.25”, (1) - 6.5”, and (1) - 9.75”
- Cart net weight: 113 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)