

Harloff Scope Storage Cabinet with Locking Easy View Doors | SC16DRDP



Stock Photo - Options May Vary

FACTORY NEW

Harloff Scope Storage Cabinet with Locking Easy View Doors: SC16DRDP. The SC16DRDP is a state-of-the-art stationary cabinet that stores and dries up to 16 scopes.

This cabinet features two fans and a HEPA filter for drying scopes. The SC16DRDP can be purchased with or without the HEPA filtration and drying fans.

Accommodates colonoscopies, EUS Scopes, TEE probes, bronchoscopes, and cystoscopes. Note: Click the following link to view additional Scope Drying Cabinets. Click the following link to view other Case Carts & Operating Room Cabinets.

FEATURES

- Drying fans 100-240VAC
- Universal Input HEPA filter captures up to 99.97% of particles, 0.3 microns, exceeding AORN scope storage recommendations
- Holds up to 16 scopes in a secure stationary location
- A lockable tambour door helps protect and secure scopes
- Engineered to protect scopes and provide easy access

SPECIFICATIONS

- Construction: Sturdy, 18 gauge steel
- Dimensions: 93.25"H x 36.10"W x 18.75"D (236.9 cm H x 91.7 cm W x 47.6 cm D)
- Cabinet Net Weight: 236 pounds (107 kg)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)