

Harloff Standard Med Cart, DISPILL USA® Compatible | SL100DSPN

FACTORY NEW

Harloff SL100DSPN Standard Med Cart, DISPILL USA® Compatible. Harloff has expanded its medication cart line to include DISPILL® medication packaging.

In addition, a standard in Canada, DISPILL®, is now being introduced in the United States. DISPILL® packaging is a unique size and requires a cart with the flexibility in design to include rails and narcotics storage specific to these dimensions. The Harloff SL100DSPN Features a Single-column medication delivery cart, DISPILL USA® Compatible.

Additionally, a Capacity of 100 plus DISPILL® cards is featured.



Stock Photo - Options May Vary

FEATURES

- Single-column medication delivery cart
- The cart has a capacity of 100 plus DISPILL® cards
- Keylock system with a Best® brand lock on the cabinet
- Illinois lock on narcotics drawer located in top 13" drawer, designed to hold DISPILL® cards
- Painted pull-out shelf

SPECIFICATIONS

- Construction: Painted steel cabinet and aluminum drawers
- Internal Drawer Dimensions: 29.5" W x 18" D
- Drawers Configuration: 1 ea - 4" multipurpose storage drawer, one ea - 3.25" multipurpose storage drawer, one ea - 13" drawer with aluminum riveted dividers and locking narcotics box., one ea - 13" drawer with aluminum riveted dividers
- Assembled Cart Dimensions: 45.75" H x 36" W x 23.5" D
- Writing surface: 44"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)