

Harloff Standard Line 36 Bin Medication Cart | SL36BIN5



Stock Photo - Options May Vary

FACTORY NEW

Harloff SL36BIN5 Standard Line 36 Bin Medication Cart. Nurses prefer medication carts with a two-column design for their ease of use. The narcotics and liquid storage in this med cart are separated from the rest of the cart.

Furthermore, this allows for easy access as needed. In addition, the drawers in this cart enable nursing staff to remove each bin as needed.

We offer both 3" and 5" bin options. The Harloff SL36BIN5 Features: Dual column medication cart with 36 - 5" wide medication bins.

FEATURES

- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Allow the nursing staff to remove each bin as needed.
- Dual-column medication cart with 36 - 5" wide medication bins. Dividers and labels are included
- BEST® brand lock on the cabinet; Illinois lock on narcotics drawer
- Painted steel pull-out shelf

SPECIFICATIONS

- Internal Drawer Dimensions: Left Side: 17" D X 27.5" W, Right Side: 17" D x 6.25" W
- Assembled Cart Dimensions: 45.75" H x 43" W x 23.5" D
- Writing Surface: 44"
- Cart Weight: 220 lbs
- Approx. Shipping Weight: 265 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)