

# Harloff Standard Line 42 Med-Bin Cart Package | SL42BIN3



Stock Photo - Options May Vary

## FACTORY NEW

Harloff SL42BIN3 Standard Line Med-Bin Cart with 42 - 3.5" Bins Standard Package. Nurses prefer medication carts with a two-column design for their ease of use. The narcotics and liquid storage in this med cart are separated from the rest. This allows for easy access as needed.

In addition, the drawers in this cart enable nursing staff to remove each bin as needed. The Harloff SL42BIN3 Features:  
Dual column medication cart with 42 - 3.5" wide medication bins. Dividers and labels are included.

BEST® brand lock on the cabinet; Illinois lock on narcotics drawer. Painted pull-out shelf. Ball-bearing, full extension drawer slides. Plastic-replaceable top with integrated push handles.

## FEATURES

- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Allows the nursing staff to remove each individual bin as needed
- Painted steel cabinet and aluminum drawers
- Dual-column medication cart
- 42 - 3.5" wide medication bins

## SPECIFICATIONS

- Internal Drawer Dimensions: Left side: 17" D x 27.5" W, Right side: 17" D x 6.625" W
- Assembled Cart Dimensions: 45.75" H x 43" W x 23.5" D
- Writing surface: 44"
- Cart weight: 220 lbs
- Approx shipping weight: 265 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)