

Harloff Standard Line 42 Med-Bin Cart Specialty Package | SL42BIN3SP



Stock Photo - Options May Vary

FACTORY NEW

Harloff SL42BIN3SP Standard Line Med-Bin Cart Specialty Package: The cart comes with a BEST® lock on the narcotics drawer, drawer tray with adjustable dividers, eye drop tray, aluminum waste container, and locking sharps container, along with the standard 42 - 3.5" Bins Standard Package SL42BIN3. Furthermore, Nurses prefer Medication carts with a two-column design for their ease of use.

The Harloff SL42BIN3 Features: This dual-column medication cart has 42 - 3.5" wide bins, with separated narcotics and liquid storage for easy access. The painted steel cabinet and aluminum drawers feature ball-bearing, full-extension slides.

It includes dividers, labels, and a single-width narcotics drawer with an Illinois lock.

FEATURES

- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Allows the nursing staff to remove each bin as needed
- Painted steel cabinet and aluminum drawers
- Dual-column medication cart
- 42 - 3.5" wide medication bins

SPECIFICATIONS

- Internal Drawer Dimensions: Left side: 17" D x 27.5" W, Right side: 17" D x 6.625" W
- Assembled Cart Dimensions: 45.75" H x 43" W x 23.5" D
- Writing surface: 44"
- Cart weight: 220 lbs
- Approx shipping weight: 265 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)