

# Harloff Standard Line 54 Bin Medication Cart w/Key Lock | SL54BIN3



Stock Photo - Options May Vary

## FACTORY NEW

The Harloff SL54BIN3 is a dual-column medication delivery cart designed for efficient nursing workflows. It features 54 medication bins organized across two columns, with narcotics and liquid storage separated from general medication storage for easy access and regulatory compliance. The cart includes removable bins, allowing nursing staff to extract individual drawers as needed during medication distribution.

The cart is built with a painted steel cabinet and aluminum drawers, featuring ball-bearing full-extension self-closing slides for smooth operation. Security is maintained through a dual locking system: a BEST® Brand Lock on the main cabinet and an Illinois Lock on the narcotics drawer, both meeting regulatory standards.

Mobility is provided by 5-inch ball-bearing premium rubber casters with full swivel capability; two casters include toe-locking brakes for stability. The cart includes a plastic replaceable top with integrated push handles and protective corner bumpers to protect both the equipment and facility infrastructure.

## FEATURES

- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Drawers in this cart allow nursing staff to remove each bin as needed
- Dual Locking System complies with regulatory security standards
- Ball-bearing, full-extension, self-closing drawer slides
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe-locking brakes, and one directional

## SPECIFICATIONS

- Assembled Dimensions: 45.75" H x 43" W x 23.5" D
- Cart net weight: 220 lbs
- Dual Column Medication Delivery Cart with 54 - 3.5" wide Medication Bins
- Internal Drawer Dimensions: Left Side: 17" D x 27.5" W, and Right Side: 17" D x 6.25" W
- Approx shipping weight: 265 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*