

Harloff Punch Card Medication Cart | SL550PC



Stock Photo - Options May Vary

FACTORY NEW

Harloff SL550PC Punch Card Medication Cart. This is a Standard Punch Card Cart with a Lock on the Cabinet and a capacity of 550 Cards. The standard package includes a pull-out shelf, adjustable punch card row distributors, and a single-column cart with an internal narcotics box.

See the tabs for more information. Are you looking for a cart with more or fewer drawers, MR-Conditional, or a cart for a specific function? Follow the links below to view a wide variety of Medical Carts, Case Carts & Operating Room Cabinets, Specialty Carts, and Popular Cart Accessories.

Have other questions? We have answers!

FEATURES

- Single-column medication cart with up to 550 punch card capacity, single-width internal narcotics box
- BEST® brand lock with four keys on the cabinet
- Illinois lock with four keys on narcotics drawer
- Painted pull-out shelf
- Ball-bearing, full-extension drawer slides

SPECIFICATIONS

- One each - 4" multipurpose storage drawer
- One each - 9.75" punch card drawer with dividers and locking narcotics box
- Two each - 9.75" punch card drawers with dividers
- 17" D x 29.5" W
- 45.75" H x 36" W x 23.5" D (116.2 cm x 91.4 cm x 59.7 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)