

# Harloff Electronic Punch Card Medication Cart | SL600PC-EKC



Stock Photo - Options May Vary

## FACTORY NEW

The Harloff SL600PC-EKC is a dual-column medication cart designed to hold 600 punch cards or blister packs with separate narcotics and liquid storage for convenient access. The cart features a CompX® electronic lock on the cabinet with BEST® brand key lock override, plus an Illinois lock on the narcotics drawer. Access rights can be assigned for up to 250 users per cart, with drawers automatically re-locking after a pre-set or programmed time interval.

The cart's two-column design improves workflow efficiency for nursing staff. It includes a painted pull-out shelf, ball-bearing full-extension drawer slides, and a plastic replaceable top with integrated push handles. The unit is equipped with 5" ball-bearing premium rubber casters—two with toe-locking brakes and one directional—for smooth mobility and stability.

Constructed with a painted steel cabinet and aluminum drawers, the cart features a hard-baked powder coat finish available in 13 colors. Protective corner bumpers provide additional durability. The cart weighs 185 lbs. and measures 45.75" H x 43" W x 23.5" D.

## FEATURES

- Dual-column medication cart with 600 punch card capacity, single-width narcotics drawer
- CompX® electronic lock on the cabinet with BEST® brand key lock override; Illinois lock on narcotics drawer
- Access rights may be assigned for up to 250 users per cart
- Drawers automatically re-lock after access (pre-set or programmed time)
- Ball-bearing, full-extension drawer slides

## SPECIFICATIONS

- Assembled cart dimensions: 45.75" H x 43" W x 23.5" D
- Writing surface: 44"
- Cart weight: 185 lbs.
- Approx. shipping weight: 230 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)