

# Harloff Standard Line Multi-Dose Medication Cart | SL72MD



Stock Photo - Options May Vary

## FACTORY NEW

Harloff SL72MD Standard Line Multi-Dose Medication Cart. Harloff multi-dose medication carts were designed to accommodate plastic bins used to hold automated strip packaging. These carts include internal drawer dividers, a pull-out shelf, and a locking narcotics box as standard features.

Nurses prefer medication carts with a two-column design for their ease of use. The Harloff SL72MD Features: Dual column unit-dose medication cart. BEST® brand lock on the cabinet, Illinois lock on narcotics drawer.

Internal drawer dividers designed to hold multi-dose bins. Painted steel pull-out shelf. Ball-bearing, full extension drawer slides. Plastic-replaceable top with integrated push handles.

## FEATURES

- Dual column unit-dose medication cart
- BEST® brand lock on the cabinet
- Narcotics and liquid storage in this med cart is separated from the rest of the cart
- Illinois lock on narcotics drawer
- Internal drawer dividers designed to hold multi-dose bins

## SPECIFICATIONS

- Internal Drawer Dimensions: Left Side: 17" D x 27.5" W, Right Side: 17" D x 6.25" W
- Assembled Cart Dimensions: 45.75" H x 43" W x 23.5" D
- Drawer Configuration: Left Column Drawers: 6 ea - 4" drawers with internal dividers included one ea - 8" drawer with two adjustable dividers. Right Column Drawers: 2 ea - 4" Multipurpose Storage Drawers, 1 - 8" Locking Narcotics Drawer, and two ea - 8" Multipurpose Storage Drawers
- Writing surface: 44" H
- Cart weight: 210 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*