

Harloff Stainless Steel 8 Drawer Cast Cart Standard Pkg | 6020



Stock Photo - Options May Vary

FACTORY NEW

Harloff Stainless Steel 8 Drawer Cast Cart Standard Pkg: 6020. This orthopedic cast cart system stores instruments and supplies for treating strains, sprains, and fractures in a cost-effective way. The cart includes a stainless steel cabinet, drawers, a bucket holder with a vinyl mat, and a 24" ruler.

Additional 6020 Features: Aluminum towel bar, Internal drawer dividers, Protective Corner Bumpers", 5" ball-bearing casters, all full swivel - one caster tracking, one with brake, and Ball-bearing full extension drawer slides. The cart comes fully assembled, is latex-free, and includes a twelve-year warranty.

Designers worked with orthopedic technicians to maximize functionality. Note: The lock in the image above is not an exact representation.

FEATURES

- Stainless steel bucket holder with a vinyl mat
- 24" stainless steel ruler mounted on top of the cart
- Aluminum towel bar
- Internal drawer dividers
- Protective Corner Bumpers

SPECIFICATIONS

- Assembled cart dimensions: 43" H x 44" W x 21" D
- Cart weight: 182 lbs.
- Approximate shipping weight: 222 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)