

# Harloff Start Cart for IV and Epidural Procedures w/Key Lock | 3123K



Stock Photo - Options May Vary

## FACTORY NEW

Harloff Start Cart for IV and Epidural Procedures w/Key Lock: 3123K. Designed to Reduce the Risk of Accidental Needle Sticks to patients and staff alike while improving the efficiency of the IV and epidural insertion procedures. Dual work surfaces and an articulating arm provide a work area to prepare and dispose of IV and Epidural supplies. This Cart also offers compact storage for point-of-use items.

Three Drawer Procedure Cart with Articulating Arm, Adjustable Push Handle, and Key Lock. Functional procedure carts improve staff efficiency and patient satisfaction. Harloff's Start Cart allows I. V.

Starts and Epidurals to be organized for administering and recovering needles, catheters, and waste.

## FEATURES

- Articulating arm allows the clinician to work from a sitting or standing position
- Dual-lipped working surfaces with a clear plastic cover overlay
- Adjustable push handle in 1-inch increments
- Soft grip handles for moving the cart from the sitting position
- 3-gallon sharps and waste containers make it easy to dispose of supplies

## SPECIFICATIONS

- Construction: Steel Cabinet and Drawers
- Base cart: 21.50"W x 22.50"D x 30" H
- Articulating arm extends height to 45.75" H
- Cart with waste and sharps containers: 31.50" W x 22.50" D
- Adjustable Push Handle: Lowest height 34.75" extends to 39.75" in 1" increments

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)