

Harloff V-Series 4 Drawer Infection Control Cart w/Key Lock | V30-4K



Stock Photo - Options May Vary

FACTORY NEW

Harloff V-Series 4 Drawer Infection Control Cart w/Key Lock. Harloff's V-Series Carts are value-priced for use in a variety of hospital departments.

This drawer configuration is often popular for treatment, procedure, and isolation applications. The cart has a Key lock with two keys, an Aluminum side rail/push handle, and 30" of vertical drawer space.

Features Include Bumper frame with projected front casters and (4) corner bumpers, 5" ball-bearing casters, all full swivel - two with brakes, Pre-drilled threaded holes for V-Series side-mounted accessories, and Ball-bearing full extension drawer slides. Additional features include Construction: Steel body and drawers, Latex-free, and a durable, hard-baked powder coat finish.

FEATURES

- Key lock with two keys
- Aluminum side rail/push handle
- 30" of vertical drawer space
- Bumper frame with projected front casters and (4) corner bumpers
- 5" ball-bearing casters, all full swivel - two with brake

SPECIFICATIONS

- Overall Cart Dimensions: 40" H x 33.25" W x 22" D
- Cart weight: 105 pounds
- Approximate shipping weight: 125 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)