

Harloff V-Series 4 Drawer ER Crash Cart w/Breakaway Lock | V24-4B



Stock Photo - Options May Vary

FACTORY NEW

Harloff V24-4B V-Series Four Drawer Emergency Crash Cart with Breakaway Lock. This crash cart is designed for healthcare facilities under budget pressure requiring a durable emergency cart offering Harloff's quality construction at a value price. Furthermore, this cart is perfect for hospitals, clinics, dialysis centers, doctors' offices, and urgent care facilities to store, transport, and dispense critical supplies used for emergencies.

The Harloff V24-4B features a break-away lock bar and a pack of 100 break-away seals. It offers 24" of vertical drawer space and an aluminum side rail/push handle.

Additionally, it includes a bumper frame with projected front casters and four corner bumpers.

FEATURES

- Break-away lock bar
- Pack of 100 break-away seals
- 24" of vertical drawer space
- Aluminum side rail/push handle
- Bumper frame with projected front casters

SPECIFICATIONS

- Internal Drawer Dimensions: 23" W x 16.75" D
- Construction: Steel body and drawers Overall Cart
- Dimensions: 34" H x 29.5" W x 22" D
- Drawer Configuration: Two-3", One-6", and One-12"
- Cart Weight: 98 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)