

# Harloff V-Series 4 Drawer ER Crash Cart w/Emergency Accessory Pkg | V24-4EMG



Stock Photo - Options May Vary

## FACTORY NEW

Harloff V24-4EMG V-Series 4 Drawer ER Crash Cart w/Emergency Accessory Pkg: This crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price. Additionally, the Crash Cart includes all commonly required accessories. This cart is latex-free and Durable.

Furthermore, it includes a twelve-year warranty, ships with wheels installed, and requires minimal assembly upon delivery. The Harloff V24-4EMG features a break-away lock bar, a pack of 100 break-away seals, and 24" of vertical drawer space.

It features a bumper frame, projected front casters, four corner bumpers, and 5" full swivel ball-bearing casters, two with brakes.

## FEATURES

- Break-away lock bar
- Pack of 100 break-away seals
- 24" of vertical drawer space
- Bumper frame with projected front casters
- (4) corner bumpers, 5" ball-bearing casters - all full swivel - two with brake

## SPECIFICATIONS

- Internal Drawer Dimensions: 23" W x 16.75" D
- Drawer Configurations: Two-3", One-6", and One-12"
- Construction: Steel body and drawers Overall Cart
- Dimensions: 36.25" H x 33.25" W x 22" D
- Height includes IV Pole

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)