

Harloff V-Series 6 Drawer ER Crash Cart w/ER Accessory Pkg | V30-6EMG



Stock Photo - Options May Vary

FACTORY NEW

Harloff V30-6EMG V-Series 6 Drawer ER Crash Cart w/ER Accessory Pkg: This crash cart is designed for healthcare facilities under budget pressure requiring a durable emergency cart offering Harloff's quality construction at a value price. Additionally, the cart includes all commonly required accessories.

The Harloff V30-6EMG Features: This cart is perfect for hospitals, clinics, dialysis centers, doctors' offices, urgent care, and nursing homes needing emergency supply carts. Furthermore, It features a break-away lock bar, 100 break-away seals, 24" vertical drawer space, and an aluminum side rail/push handle.

The bumper frame has projected front casters, four corner bumpers, and 5" ball-bearing casters, all full swivel—two with brakes.

FEATURES

- Break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Painted steel top rail/push handle
- Bumper frame with projected front casters and (4) corner bumpers

SPECIFICATIONS

- Internal Drawer Dimensions: 23" W x 16.75" D
- Construction: Steel Body and Drawers
- Overall Cart Dimensions: 67.4" H x 33.25" W x 22" D
- Height includes IV Pole
- Cart weight: 138 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)