

Harloff Value Line Punch Card Medication Cart w/Key Lock | VLT360EPC



Stock Photo - Options May Vary

FACTORY NEW

Harloff VLT360EPC Value Line Punch Card Medication Cart with Key Lock. The Harloff Value Line medication cart series engineers a cart that meets the needs of smaller facilities and nursing units.

This economically priced cart includes a pull-out shelf and a locking narcotics box as standard features. Additionally, the cart design accommodates punch cards or blister packs.

Are you looking for more/fewer drawers, MR-Conditional, or a cart for a specific function? Follow the links below to view a wide variety of Medical Carts, Case Carts & Operating Room Cabinets, and Specialty Carts & Accessories. Have other questions? We have answers!

FEATURES

- Cart style is designed for punch cards or blister packs
- Tall cabinet with key locking system
- Taller drawer height for increased capacity
- Has a capacity of up to 360 punch card carts
- Internal locking narcotics box

SPECIFICATIONS

- Construction: Painted Steel Cabinet and Drawers.
- Assembled Cart Dimensions: 46.5" H x 29.5" W x 22.5" D
- Writing surface: 42.5"
- Drawer configuration: 1 - ea 3" multipurpose storage drawer, 2 - ea 9.75" drawers with two adjustable dividers, and 1 - ea 9.75" drawer with one divider and locking narcotics box.
- Cart weight: 140 pounds

**Brochure and images are a general overview of the product. Actual product options may vary.*