

Harloff Wood Vinyl 4 Drawer Treatment Cart, Standard Package | WV6350



Stock Photo - Options May Vary

FACTORY NEW

Harloff Wood Vinyl 4 Drawer Treatment Cart, Standard Package: WV6350-CM. Wood-style treatment carts provide the same storage as standard carts, with a more welcoming and relaxing appearance.

The WV6350 is a four-drawer treatment cart made from lightweight aluminum with a wood vinyl coating, designed to resemble furniture. Each cart features a removable top with a granite look and can come in any of our wood vinyl colors.

This model is shown in our standard Cherry Mahogany color.

FEATURES

- This cart is designed to look like a piece of furniture Lightweight vinyl-coated aluminum
- Key lock with two keys
- 5" (12.7cm) ball-bearing premium casters, all full swivel, one caster braking, one-directional
- Ball-bearing full extension drawer slides
- Pull out shelf

SPECIFICATIONS

- Vinyl-coated aluminum
- The top is MDF covered with a heat-pressed thermal foil laminate
- two ea - 3.25" drawers
- one ea - 6.5" drawer
- one ea - 12" drawer

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)