

# Harloff Wood Vinyl 240 Capacity Punch Card Medication Cart, Std Package | WV240PC



Stock Photo - Options May Vary

## FACTORY NEW

Harloff WV240PC Wood Vinyl 240 Capacity Punch Card Medication Cart with Standard Package. This wood blend medication cart holds up to 240 punch cards and supports Culture Change initiatives in long-term care by providing a home-like feel. Constructed from lightweight aluminum with a wood vinyl coating, this four-drawer cart features a removable top with a granite look.

Additional Information: The cart comes in any of our wood finishes and is shown in our standard Cherry Mahogany. Each wood blend medication cart features a standard key lock, adjustable punch card row dividers, ball-bearing full extension drawer slides, protective corner bumpers, and 5" ball-bearing premium rubber casters.

## FEATURES

- Lightweight vinyl-coated aluminum
- Single-column medication cart with 240 punch card capacity
- Key locking cabinet with internal locking narcotics box located in the top drawer
- Adjustable row dividers for punch cards and liquid storage
- Ball-bearing, full-extension drawer slides

## SPECIFICATIONS

- one each - 3" multipurpose storage drawer
- one each - 9.75" punch card drawer with locking narcotics and adjustable dividers
- two each - 9.75" punch card drawer with adjustable dividers
- 44.5" H x 19" W x 25" D
- Writing Surface: 44"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)