

# Harloff Wood Vinyl 540 Capacity Punch Card Med Cart Pkg | WV540PC



Stock Photo - Options May Vary

## FACTORY NEW

Harloff WV540PC Wood Vinyl 540 Capacity Punch Card Medication Cart, Std Package. Culture change initiatives have increased the interest in wood medication carts for use in long-term care. These mobile furniture-style med carts offer the look of furniture along with the function, security, and narcotics storage that facilities require.

The Harloff WV540PC Wood Look Med Cart offers various sizes and drawer options for different medication packaging. Made from lightweight, vinyl-coated aluminum, it holds 540 punch cards with a key-locking cabinet and narcotics box in the top drawer.

The cart features adjustable dividers, full-extension ball-bearing drawer slides, a pull-out shelf, and 5" premium rubber casters.

## FEATURES

- Vinyl-coated aluminum
- Lightweight vinyl-coated aluminum
- Single-column medication cart with 540 punch card capacity
- This cart is designed to look like a piece of furniture
- Key locking cabinet with narcotics box located in the top drawer

## SPECIFICATIONS

- Assembled Cart Dimensions: 45" H x 34.5" W x 24.5" D
- Writing Surface: 44"
- Cart Weight: 127 lbs
- Approx. Shipping Weight: 157 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)