

# Hassoc Ergonomic Lab / Anesthesia Chair | 8302



Stock Photo - Options May Vary

## FACTORY NEW

Hassoc Ergonomic Lab / Anesthesia Chair: 8302. Hassoc's Ergonomic Lab Chair features an Aluminum Base with/ Toe Caps, an Adjustable Seat Back, an 18" diameter footrest standard with casters shown, and the 3/4" thick commercial grade seat virtually eliminates warping and breakage.

No bottoming out on the seat pan! Soft Wheel Casters or Glides are Optional. The 8302 comes with a Limited Lifetime Warranty.

## FEATURES

- Aluminum Base w/ Toe Caps
- Adjustable Seat Height
- 18" Diameter Footrest Standard with Casters
- Shown Soft Wheel Casters or Glides Optional
- A 3/4" thick commercial-grade seat virtually eliminates warping and breakage

## SPECIFICATIONS

- Base Construction: Aluminum Base
- Seat Height Adjustment: 23" to 33"
- Seat Width: 20"
- Seat Cushion: 3.5" thick foam
- Base Diameter: 24"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)