

# Hassoc Hydraulic Surgical Stool | 7404



Stock Photo - Options May Vary

## FACTORY NEW

Hassoc Hydraulic Surgical Stool: 7404. Hydraulic Surgical Stool Seat Height Adjustment: 19.5" - 26.5" Seat Diameter: 15.75" High. This Surgical Stool has Density Fire-Resistant Foam and 3.5" thick foam. It provides a firm but comfortable seat with no bottoming out on the seat pan! The 3/4" thick commercial grade seat virtually eliminates warping and breakage.

Hassoc 7404 Additional Features: Chrome 1/2" Solid Steel Foot Control. Rotation Lock on the seat. The Use of ISO 9002 components for most metal parts guarantees quality. Waterfall Seam for Cleaning Ease. Black Aluminum Base.

Limited Lifetime Warranty. Available Options: Back Rest. Black or Silver Powder Coat Hydraulic Finish.

## FEATURES

- Seat Height Adjustable
- Density Fire-Resistant Foam
- 3.5" Foam Seat
- Solid Steel Foot Control
- Seat Rotation Lock

## SPECIFICATIONS

- Seat Height Adjustment: 19.5" - 26.5"
- Seat Diameter: 15.75"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)