

Hassoc Hydraulic Surgical Stool | 7408



Stock Photo - Options May Vary

FACTORY NEW

The Hassoc Hydraulic Surgical Stool Model 7408 features a hydraulic seat height adjustment system designed for surgical and clinical environments. The stool incorporates a 3.5-inch thick density fire-resistant foam seat that provides firm, comfortable support without bottoming out, paired with a 3/4-inch thick commercial-grade seat pan that resists warping and breakage.

The stool is equipped with a chrome 1/2-inch solid steel foot control for smooth height adjustment and includes a rotation lock to stabilize the seat position. The brushed aluminum base provides durability and easy mobility in clinical settings.

Construction utilizes ISO 9002 components for metal parts to ensure quality standards. The waterfall seam design facilitates cleaning and maintenance. The stool comes with a limited lifetime warranty and is available with optional back rest and black or silver powder coat hydraulic finish.

FEATURES

- Seat Height Adjustable
- Density Fire-Resistant Foam
- 3.5" Foam Seat
- Solid Steel Foot Control
- Rotation Lock

SPECIFICATIONS

- Seat Height Adjustment: 19.5" - 26.5"
- Seat Diameter: 15.75" H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)