

Hausted 4160 Pediatric Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Hausted 4160 Pediatric Stretcher provides 360° safety while maintaining complete patient access. Hydraulically operated height adjustment ranges from 32.25" to 39", with dual side foot pedals for easy operation. Trendelenburg capability and a manually adjustable backrest accommodate various clinical positioning needs.

The stretcher features removable end rails for unrestricted patient access and rolls smoothly on 8" casters. Brake and steering functions can be activated from all four corners of the base, providing convenient control from any position. The polymer-coated full-length vertical side rails and full perimeter non-marring bumper system protect both patient and equipment.

Overall dimensions are 30" W x 60" L with a patient surface of 25.25" W x 56" L. The plastic base includes built-in oxygen tank holder and storage. Weight capacity is 325 lbs, with IV pole receptacles and accessory wells at each corner.

FEATURES

- Dual Side Foot Pedal Hydraulic Controls
- Trendelenburg
- Manual Backrest Adjusts with One Hand
- Brake and Steering Activation from All 4 Corners
- Plastic Base with Storage and Built-In Oxygen Tank Holder
- Full Perimeter Non-Marring Bumper System
- Polymer Coated Removable End Rails Provide Unrestricted Patient Access
- Weight Capacity: 325 lbs

SPECIFICATIONS

- Height Range: 32.25" to 39"
- Trendelenburg: 0° - 10
- Manual Adjustable Backrest: 0° - 45° with One Hand
- Patient Surface: 25" x 56"
- Side Rails: 18" x 55"
- Weight Capacity: 325 lbs
- Bumper System: Full Perimeter Non-Marring
- IV Pole Receptacles: 1 at Each Corner

**Brochure and images are a general overview of the product. Actual product options may vary.*