

Pedigo 500 Pediatric Crib – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Pedigo 500 Pediatric Crib is a completely refurbished crib/stretcher featuring a unique 2-post oval design that rolls easily on 5" locking casters. The crib offers one-handed operation for raising and lowering side rails with an intermediate stop position, and features 2" bar spacing for safety. The patient surface measures 26" wide by 56" long.

Key positioning capabilities include three-position locking Fowler, Trendelenburg/Reverse Trendelenburg adjustment, and easy side rail adjustment. The counterbalanced side rails prevent sudden drops and allow caregivers to access patients safely while holding children. Two 5" locking steer casters with central locks provide maximum maneuverability.

This model remains available as a current production unit and functions as both a pediatric crib and stretcher. A new mattress is included. Optional accessories such as a basic IV pole, croup top, and utility shelf are available for separate purchase.

FEATURES

- Completely, Professionally Refurbished
- One-Handed Side Rail: Raise and Lower
- Side rail intermediate stop position
- 2" Rail Spacing
- Trendelenburg / Reverse Trendelenburg
- Full Function Fowler
- 2 Central Locking, 5" Casters
- six months parts Warranty

SPECIFICATIONS

- Top Rail Height: 59 1/2"
- Height with Top: 74"
- Length: 60"
- Width: 30 1/4"
- Spring Height: 32 1/4"
- Caster Diameter: 5"
- Fowler Angles: 0,15,30,45 degrees
- Trendelenburg Angles: -5,0,5,10 degrees

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)