

Replacement Side Rail Bumper Pads for Hausted 4160 Pediatric Crib



Stock Photo - Options May Vary

FACTORY NEW

Replacement side rail bumper pad set for the Hausted 4160 Pediatric Crib. This set includes two protective pads designed to cushion and protect the crib's side rails. The pads are constructed with high-density foam and covered in durable Herculite Sure Check Comfort Fabric.

The pads feature welded corners to prevent fluid seepage and are treated for flame and fluid resistance to meet safety requirements. The fabric is resistant to bacteria, stains, and tears, ensuring long-term durability in clinical environments.

Each pad measures 1 inch by 15.5 inches by 55 inches and comes in black. End bumper pads are sold separately.

FEATURES

- Side Rail Bumper Pads: Set of two
- Pad thickness: 3 inches
- Color available: Black
- High Density foam
- Covered in durable Herculite Sure Check Comfort Fabric
- Welded corners to prevent fluid seepage
- Flame and Fluid resistance for safety
- Resistance to bacteria, stains and tears

SPECIFICATIONS

- (2) 1" x 15.5" x 55" Pads contain High Density foam

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)