

# Stryker Cub FL19H / 190 Pediatric Crib - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker Cub Pediatric Crib (Model FL19H / 190) is designed to support safe and effective care for pediatric patients. It features precise hydraulic adjustments for height, Trendelenburg, and reverse Trendelenburg positioning, along with central brake and steer functions that operate foot-actuated from the sides for hands-free operation. The crib includes a manually adjustable backrest with integrated foot prop support.

The crib offers counterbalanced siderails for easy operation and unobstructed access when lowered. Additional features include swing open access doors, a built-in oxygen tank holder, a retractable 5th wheel for maneuverability, and 6" casters. The strong build and child-friendly design make it suitable for pediatric care environments.

This refurbished unit has been thoroughly inspected for quality and comes with a 6-month warranty.

## FEATURES

- Retractable Backrest
- Swing Open Doors
- Easily-Lowerable Side Rails for full access
- Manual Foot Pump
- O2 Tank Holder
- Central Brake and Steer
- Hands-Free Trendelenburg/Reverse Trendelenburg
- Warranty

## SPECIFICATIONS

- Overall Length: 65.5"
- Overall Width: 37"
- Height Range: 33" - 40"
- Weight Capacity: 400 lbs
- Backrest Angle: 0° / 25° / 40° / 55°
- Trend/Rev Trend: +/- 9.5°
- Patient Surface: 29.5" x 57"
- Includes New 3" Pad

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)