

Medtronic Physio-Control LIFEPAK 20 Defibrillator - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Medtronic Physio-Control LIFEPAK 20 is a defibrillator/monitor designed for hospitals and clinics, combining 3-lead biphasic defibrillation, pacing, and AED capabilities. The device features an intuitive interface with the Shock Advisory System that guides users through voice and visual prompts in a simple 1-2-3 operation, making it accessible for infrequent AED-trained responders.

With an easy-push latch, the LIFEPAK 20 converts from AED mode to manual defibrillator mode, displaying waveforms and providing access to advanced cardiovascular life support tools. The device supports non-invasive pacing, ECG monitoring with 3- or 5-wire configurations, synchronized cardioversion, and optional Masimo SET pulse oximetry (SpO2).

Compact and lightweight, the LIFEPAK 20 is designed for indoor use and can be quickly deployed to the scene or used during patient transport.

FEATURES

- 3-lead, Biphasic, Pacing, and AED
- Easy To Use - Ideal for infrequent AED-trained responders to quickly understand and use
- "Shock Advisory System™" guides the user with voice and visual prompts through each step and simple 1-2-3 operation
- Automatically converts to a manual defibrillator, including waveform displays
- In manual mode, all of the advanced cardiovascular life support tools become available
- The 20 offers non-invasive pacing, Masimo SET pulse oximetry-SpO2 (optional), ECG monitoring (3- or 5-wire), and synchronized cardioversion
- Compact, lightweight, and easy to rush to the scene or use during transport
- Created especially for hospitals and clinics

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com