

Zoll M-Series Defibrillator Biphasic | Pacing | AED - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Zoll M-Series Defibrillator Biphasic | Pacing | AED: Refurbished. The Zoll M-Series defibrillator is one of the smallest and most lightweight defibrillators available in its class. This unit features a biphasic waveform, non-invasive pacing, and an AED advisory mode.

Operating the system is simple and user-friendly, with straightforward controls and easy-to-use menus. Having doubts about what is the correct procedure to check your Defib every day? Check out our blog post: '6 Tips on Checking Your Defibrillator Daily,' to learn more about it.

To view additional units and accessories, click the following link to browse Venture Medical's Complete Catalog of New and completely refurbished defibrillators and accessories. This Product is Discontinued.

FEATURES

- Biphasic Waveform Technology
- Pacing Function
- AED Function
- Synchronized Cardioversion
- 3 Lead EKG

SPECIFICATIONS

- Max Joules: 360j
- Biphasic up to 200j
- Charge Time: Less than 7 seconds
- Display Size (Screen): 14.2 cm (5.6")
- Monitor Displays Selected and Delivered Energy

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)