

Medtronic Physio-Control LIFEPAK 20e Defibrillator - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Medtronic Physio-Control LIFEPAK 20e is a defibrillator/monitor designed for hospital and clinic crash carts. It features 3-lead biphasic defibrillation, pacing, and AED capabilities with an intuitive interface suitable for infrequent AED-trained responders. The device includes a battery status indicator on the display and enhanced battery management technology.

The LIFEPAK 20e guides users through each step with voice and visual prompts via the Shock Advisory System™, enabling simple 1-2-3 operation. It automatically converts to manual defibrillator mode with waveform displays, making all advanced cardiovascular life support tools available to trained operators.

The unit offers non-invasive pacing, optional pulse oximetry (SpO2), ECG monitoring with 3- or 5-wire configurations, and synchronized cardioversion. Its compact, lightweight design makes it suitable for transport and emergency response situations.

FEATURES

- 3-lead, Biphasic, Pacing, and AED
- Easy To Use - Ideal for infrequent AED-trained responders to quickly understand and use
- "Shock Advisory System™" guides the user with voice and visual prompts through each step and simple 1-2-3 operation
- Automatically converts to a manual defibrillator, including waveform displays
- In manual mode, all of the advanced cardiovascular life support tools become available
- The 20e offers non-invasive pacing, pulse oximetry-SpO2 (optional), ECG monitoring (3- or 5-wire), and synchronized cardioversion
- Compact, lightweight, and easy to rush to the scene or use during transport
- The LIFEPAK 20e defibrillator/monitor is the ideal crash cart defibrillator

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com