

Zoll R Series Plus Defibrillator – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Zoll R Series Plus is a hospital-grade defibrillator and patient monitor designed for crash carts and critical care environments. It provides fast access to essential patient data and therapy controls with an intuitive interface and automated readiness checks to help medical staff respond quickly during sudden cardiac arrest emergencies.

Advanced technologies like Real CPR Help® and See-Thru CPR® guide responders to deliver effective compressions with minimal CPR interruptions. Real-time feedback helps maintain proper compression depth and rate, improving resuscitation quality. The system powers up in AED mode for rapid deployment and supports both adult and pediatric patients with tailored CPR analysis and therapy options.

The R Series Plus offers comprehensive monitoring and device management capabilities. OneStep™ electrodes simplify defibrillation, pacing, and ECG monitoring using a single cable connection. Automated self-testing checks over 100 system indicators to ensure code-readiness. Optional monitoring for ECG, SpO₂, EtCO₂, and NIBP provides clinicians with critical patient data for effective care.

FEATURES

- Real CPR Help® provides real-time compression feedback
- See-Thru CPR® reduces interruptions during rhythm analysis
- OneStep™ electrodes for defibrillation, pacing, and ECG monitoring
- Code-Ready® automated daily self-testing
- AED and manual defibrillation modes
- Transcutaneous pacing capability
- Optional SpO₂, EtCO₂, and NIBP monitoring
- Event data recording and storage

SPECIFICATIONS

- Energy Selections: 1 to 10, 15, 20, 30, 50, 75, 100, 120, 150, and 200 joules selected using controls on the front of the defibrillator or sternum paddles. (Note: When using appropriate pediatric resuscitation electrodes, the 75-joule setting is replaced by 70- and 85-joule settings.)
- Charge Time: Less than 7 seconds to maximum energy
- Display Size (Screen): 6.5" (16.5 cm) diagonally color LCD
- Battery: 10.8 V rechargeable lithium-ion battery
- Size: 8.2" H x 10.5" W x 12.5" D
- Weight (Unit Only): 13.6 lbs (6.2 kg)
- ECG Monitoring: 3-lead or 5-lead ECG monitoring capability
- Operating Temperature: 0°C to 40°C (32°F to 104°F)

Toll-Free Local
 800-627-3215 | 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com