

# Cardiac Science Powerheart G3 AED - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Cardiac Science Powerheart G3 is a biphasic automated external defibrillator designed for rapid response in cardiac emergencies. The device features one-button operation and rescue-ready technology, allowing clinicians to quickly confirm the unit is prepared for use. Its lightweight and compact design makes it suitable for various clinical and public access settings.

The Powerheart G3 includes Rescue Coach voice prompts that guide users through the resuscitation process with clear, authoritative instructions. The voice guidance remains audible in noisy environments, helping caregivers maintain focus and confidence during critical moments.

Each refurbished unit arrives completely patient-ready and includes new adult therapy pads, a new lithium-ion battery, and a user manual. Semi-automatic or fully automatic operation modes are available to meet specific facility requirements.

## FEATURES

- Semi or Fully Automatic Units Available
- Biphasic
- Lightweight
- Compact
- New Adult Pads (1 Set)
- New Battery
- User Manual
- Completely Refurbished and Patient Ready

## SPECIFICATIONS

- Height: 3.3"
- Width: 10.6"
- Depth: 12.4"
- Weight: 6.6 lbs
- Energy Range: 150J -360J, Escalating
- Waveform: Rectilinear Biphasic
- Power: Lithium-Ion Battery
- CPR: Metronome for Compression Frequency

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)