

HeartSine Samaritan PAD AED - NEW



Stock Photo - Options May Vary

FACTORY NEW

The HeartSine Samaritan PAD AED is designed for ease of use by rescuers of all experience levels. A graphical interface with three simple action icons and clear voice instructions guides users through each step of the rescue process. The device weighs only 2.4 lbs for convenient transport.

Operation is straightforward with an auto self-test function and illuminated icons indicating "Attach pads," "Do not touch," and "Safe to touch." An integrated metronome coaches CPR compression rates during resuscitation efforts. The device uses SCOPE Biphasic waveform technology that adjusts to patient impedance for optimized shock delivery.

The unit records all events including ECG traces, shock data, and CPR intervals with time and date stamps. Data can be downloaded to a PC using optional SAVER EVO software. The Samaritan PAD is suitable for home, business, and public settings and includes a carrying case, pad-pak cartridge with battery, data port, quick instruction card, user manual, and poster.

FEATURES

- Simple operation
- A graphical interface that "illuminates" to guide the user
- Clear voice instructions
- Auto self-test
- Integrated Metronome prompts and coaches CPR chest compression rate
- Events are recorded, with the time, date, ECG trace & duration, shock delivery information, and CPR intervals
- Only 1.1 kg (2.4 lbs)
- Sophisticated SCOPE™ Biphasic waveform technology adjusts patient impedance to maximize effective shock delivery

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)