

# Hill Rom P8050 OB-GYN Stretcher – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Hill-Rom P8050 is an OB-GYN stretcher designed for fast, easy, and comfortable patient positioning during exams and procedures. It features integrated foot supports and a sliding top that allows for quick adjustments and improved access for healthcare providers. The foot section easily retracts to reveal footrests, central brake and steer system, push handles, and a built-in IV pole.

The stretcher's thoughtful design improves workflow and provides medical staff quick access to essential features needed for safe and efficient care. Gas spring assist backrest, hydraulic height operation, and Trendelenburg positioning capabilities support a variety of procedural needs.

The P8050 is suitable for clinics, surgical centers, and hospitals. It delivers long-lasting performance and patient comfort in a cost-effective package, handling many types of OB-GYN procedures while maintaining quality and reliability.

## FEATURES

- Trendelenburg/Reverse Trendelenburg
- Back Rest Up/Down (Gas Spring Assist)
- Height Up/Down
- Collapsible Side-rails
- Central Brake and Steer
- Steering
- O2 Tank Well
- Includes New Pad Set

## SPECIFICATIONS

- Overall Length / Width: 80.5" / 35.5"
- Minimum / Maximum Bed Height: 24" / 38"
- Side Rail Length: 37"
- Caster Size: 8"
- Weight Capacity: 700 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)