

# Stryker 1061 Gynnie OB/GYN Stretcher - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The Stryker 1061 Gynnie is a specialized OB/GYN stretcher designed for use in Emergency Departments, Labor and Delivery units, and gynecological procedure rooms. Fully refurbished to meet high-performance and safety standards, it provides reliable support for both patients and clinicians.

Key features include a Pneumatic Assist Fowler backrest for easy positioning, a folding leg section with integrated stirrups for improved patient access during procedures, and Trendelenburg/Reverse Trendelenburg capability. The stretcher incorporates Glideaway side rails, central brake and steering, and an integrated O2 tank well.

The stretcher offers hydraulic height adjustment ranging from 22.5" to 36", accommodates patients up to 500 lbs, and comes with a new pad and warranty.

## FEATURES

- Height: Up/Down
- Backrest: (Fowler) Pneumatic Assist
- Trendelenburg /Reverse Trendelenburg +/-
- Glideaway™ Side Rails
- Integrated O2 Tank Well
- Folding Foot Section and Stirrups
- Central Brake and Steering
- Completely Refurbished

## SPECIFICATIONS

- Overall Length: 81.5"
- Overall Width: 31.5" (Side Rails Up)
- Overall Width: 26.75" (Side Rails Down)
- Height: (Maximum) 36"
- Height: (Minimum) 22.5"
- Patient Surface: 26" x 77"
- Side Rails: 13" x 39"
- Weight Capacity: 500 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)