

Hospira Abbott Plum A+ Infusion Pump - Refurbished

FULLY REBURBISHED · CLINICALLY CERTIFIED

Hospira Abbott Plum A+ Infusion Pump - Refurbished. The refurbished Hospira Abbott Plum A+ Infusion Pump is a reliable, dual-channel infusion system designed for precise and flexible IV medication delivery. Moreover, fully refurbished and biomedically inspected, this pump meets safety standards, providing a reliable, cost-effective option for hospitals and clinics.

In addition, key features include secondary infusion, dual medication delivery, and air-in-line detection to support safe, accurate treatment. The Plum A+ system also supports a wide range of infusion therapies, including critical care, general medicine, and oncology treatments.

As a result, the Plum A+ Pump streamlines medication delivery and reduces errors with its simple interface and reliable performance.



Stock Photo - Options May Vary

FEATURES

- KVO Rate 1.0 mL/HR
- Volume Range: 0.1 to 99.9 mL (in 0.1 mL increments)
- New Rechargeable Battery
- Warranty
- Completely Refurbished

SPECIFICATIONS

- Weight: 9.5 lbs
- Dimensions: 8"x 8"x 6" (excluding pole clamp)
- Battery Life: Approx. 6 Hours
- Volume Range: 0.1 to 99.9 mL (in 0.1 mL increments)
- KVO Rate: 1.0 mL/HR or the last primary delivery rate, whichever is less

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)