

# Hospira Abbott Plum A+ Triple Infusion Pump - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Hospira Abbott Plum A+ Triple Infusion Pump - Refurbished. This Infusion Pump has been completely refurbished and bio-medically inspected.

The A+ Infusion Pump features automated second-line delivery, concurrent delivery of 6 medications, air-in-line detection, and elimination. Additional Features Include Programmable standby/delayed start settings Flow rate calculation Hospira's Clave closed system Includes New Battery, User Manual, and Warranty.

Click the following link to view Venture Medical's Complete Catalog of Infusion, Syringe, and Irrigation Pumps, Stands, and Accessories.

## FEATURES

- Delivers Infusions with Confidence with the Unique PlumSet Technology Automated Second Line Delivery
- Programmable Delayed Start
- Triple Channel (6 Concurrent medications)
- Rechargeable Battery
- Warranty

## SPECIFICATIONS

- Weight: 28 lbs with (3) Batteries
- Dimensions: 19" x 15" x 14" (including pole clamp)
- Battery Life: Approx. 6 Hours
- Mounting: Pole
- Power Requirement: 120v

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)